

Drying oven

HUMi-DRY 900

+30...+150°C



Standard equipment:

- ✓ Full AISI 304 stainless steel exterior
- ✓ AISI 304 stainless steel interior
- ✓ 2 controllers
- ✓ 4 swivel casters, 2 of which have brakes
- ✓ Adjustable pressure bulkhead for optimum drying process
- ✓ CE marking
- ✓ Mechanical espagnolette door lock
- ✓ Protection class 2.0 according to DIN 12880
- ✓ Electronic temperature controller with digital display with the possibility of changing the temperature setting every 0.1K
- ✓ ECO mode when humidity SP is achieved
- ✓ 6 stainless steel shelves

Optional equipment:

- ✓ ATEX Safety
- ✓ Power failure alarm
- ✓ Temperature sensor
- ✓ Changing the standard temperature range
- ✓ Access port
- ✓ 150 mm drawer
- ✓ Additional stainless steel shelves (shelves or stainless steel wire)
- ✓ LED lighting

HUMi-DRY 900



Uniform temperature



Electric heating



Ventilated airflow



4 swivel casters, 2 of which have brakes




Optimal drying



Customisation possible



Eco mode



Technical specifications:

	Width	Depth	Height
External dimensions [mm]	1185	838	1983
Internal dimensions [mm]	1038	603	1420
Type of heating	Electric		
Type of ventilation	Compressed		
Gross capacity	888 [liters]		
Net capacity	854 [liters]		
Weight	280 [kg]		
Max. load	300 [kg]		
Number of shelves	Max. 19		
Number of guides / horizontal supports	19		
Max. shelf load	50 [kg]		
Temperature range	+5 above ambient temp. ...+150 [°C].		
Power	400/50 [V/hz]		
Number of external doors	2		
Power rating	User-dependent		
Rated current	3500/6500/9500 [W]		